

## Custom Varicella Zoster (VZV) Inoculum

Description: VZV (Ellen) from a 6 roller bottle infection of NHDF harvested and cryopreserved to 25 x 1 ml amps.

Virus Strain: Ellen – traceable to ATCC VR-1367

Cryopreservation Media: EMEM 10% FBS 8% DMSO

Lot Number: 28-077

Quantity: 25 x 1 ml – Store in LN2

### TCID Testing using CPE or IFA (anti-gE Mab) at 7 days PI

Dilution	-1	-2	-3	-4	-5	-6	-7	-8
CPE	++++	++++	++++	++++	+++	+++	++	+
IFA	++++	++++	++++	+++	++	++	+	+/-

